



X-ray's prediction

UNIVERSITY OF PIRAEUS – DEMOCRITUS INSTITUTE

ARTIFICIAL INTELLIGENCE, MASTER

MACHINE LEARNING, COURSE

THEODOROS ADAMANTIDIS, TN2001

ATHENS, FEBRUARY 2021

Content

- Introduction
- Feature extraction
- Experiments
- Results


Introduction



- Detect pneumonia on X-Ray images.
- 5,863 X-Ray images of two categories
- Chest X-ray images (anterior-posterior) were selected from retrospective cohorts of pediatric patients of one to five years old from Guangzhou Women and Children's Medical Center, Guangzhou. All chest X-ray imaging was performed as part of patients' routine clinical care.

Feature extraction (1/3)

- Corresponding notebook **ExtractFeatures.ipynb** in GitHub.
- Developed and support 3 different types of image extraction
 - Plain flatten image array
 - KAZE descriptors array
 - HOG descriptors array
- All images resized to 100x100 size.
- Data are stored after extraction into pickle files for fast accessing in the experiments.



Feature extraction (2/3) - KAZE


- Detect and get the first 32 key points of the image ordered by their response value.
- Compute descriptor using these key points constructing a vector of 64 size.
- In case of the vector is less than 64, fill it with zeros.

Feature extraction (3/3) - HOG

Input image



Histogram of Oriented Gradients




- Cell size 10x10
- Block size 2x2
- Number of orientation bins 9
- Calculates descriptor and returns the array flatten.

Experiments (1/7)

- For the experiments used the pickle files extracted from the feature extraction step.
- Standardize features.
- Used to measure how “good” the model estimators are
 - Accuracy
 - Confusion matrix
 - Precision score
- Because of the quite difference in the number of each classes data, the accuracy by it self is not a good estimator.
- The models used are
 - Logistic regression
 - KNN
 - Discriminant analysis
 - Gaussian naive bayes
 - SVM
- Used GridSearchCV with RepeatedStratifiedKFold to tune models

Experiments (2/7)



- Number of samples in training data: 5216
 - Number of samples in training data where category is Pneumonia : 3875
 - Number of samples in training data where category is Normal: 1341
- Number of samples in testing data: 624
 - number of samples in testing data where category is Pneumonia : 390
 - number of samples in testing data where category is Normal: 234
- Number of samples in validate data: 16
 - number of samples in validate data where category is Pneumonia : 8
 - number of samples in validate data where category is Normal: 8

Experiments (3/7) – Logistic regression

No tuned

- Accuracy LR classifier on training set: 1.00
- Accuracy of LR classifier on test set: 0.80
- Accuracy of LR classifier on validate set: 0.50
- Precision score 0.7661290322580645

Pred/Label	Normal	Pneumonia
Normal	118	116
Pneumonia	10	380

Tuned (C = 0.01)

- Accuracy LR classifier on training set: 1.00
- Accuracy of LR classifier on test set: 0.80
- Accuracy of LR classifier on validate set: 0.50
- Precision score 0.768762677484787

Pred/Label	Normal	Pneumonia
Normal	120	114
Pneumonia	11	379

Experiments (4/7) – KNN

No tuned

- Accuracy of K-NN classifier on training set: 0.96
- Accuracy of K-NN classifier on test set: 0.81
- Accuracy of K-NN classifier on validate set: 0.50
- Precision score 0.8096330275229358

Pred/Label	Normal	Pneumonia
Normal	151	83
Pneumonia	37	353

Tuned (neighbors = 11, weight = 'distance', alg = kd_tree)

- Accuracy of K-NN classifier on training set: 1.00
- Accuracy of K-NN classifier on test set: 0.82
- Accuracy of K-NN classifier on validate set: 0.50
- Precision score 0.768762677484787

Pred/Label	Normal	Pneumonia
Normal	158	76
Pneumonia	37	353

Experiments (5/7) – Discriminant analysis

No tuned

- Accuracy of LDA classifier on training set: 1.00
- Accuracy of LDA classifier on test set: 0.79
- Accuracy of LDA classifier on validate set: 0.50
- Precision score 0.766260162601626

Pred/Label	Normal	Pneumonia
Normal	119	115
Pneumonia	13	377

Tuned (solver = 'lsqr')

- Accuracy of LDA classifier on training set: 1.00
- Accuracy of LDA classifier on test set: 0.79
- Accuracy of LDA classifier on validate set: 0.50
- Precision score 0.7676348547717843

Pred/Label	Normal	Pneumonia
Normal	122	112
Pneumonia	20	370

Experiments (6/7) – Gaussian naive bayes

- Accuracy of GNB classifier on training set: 0.89
- Accuracy of GNB classifier on test set: 0.76
- Accuracy of GNB classifier on validate set: 0.50
- Precision score 0.780373831775701

Pred/Label	Normal	Pneumonia
Normal	140	94
Pneumonia	56	334

Experiments (7/7) – SVM

No tuned

- Accuracy of SVM classifier on training set: 1.00
- Accuracy of SVM classifier on test set: 0.79
- Accuracy of SVM classifier on validate set: 0.50
- Precision score 0.7589641434262948

Pred/Label	Normal	Pneumonia
Normal	113	121
Pneumonia	9	381

Tuned (C = 100, kernel = 'Linear')

- Accuracy of SVM classifier on training set: 1.00
- Accuracy of SVM classifier on test set: 0.80
- Accuracy of SVM classifier on validate set: 0.50
- Precision score 0.7703252032520326

Pred/Label	Normal	Pneumonia
Normal	121	113
Pneumonia	11	379

Results

- Model:: **Logistic regression** achieved:: 0.7661290322580645 precision score.
- Model:: **Logistic regression tuned** achieved:: 0.768762677484787 precision score.
- Model:: **KNN** achieved:: 0.8096330275229358 precision score.
- Model:: **KNN tuned** achieved:: 0.8228438228438228 precision score.
- Model:: **Discriminant analysis** achieved:: 0.7676348547717843 precision score.
- Model:: **Discriminant analysis tuned** achieved:: 0.766260162601626 precision score.
- Model:: **Gaussian naive bayes** achieved:: 0.780373831775701 precision score.
- Model:: **SVM** achieved:: 0.7589641434262948 precision score.
- Model:: **SVM tuned** achieved:: 0.7703252032520326 precision score.